



NEVIFILE

An accessory to Nevisense, and a software application that you can install on your computer. NeviFile is required in order to print reports from Nevisense via the network, immediately when an examination is completed. NeviFile can also be used to save reports on a computer for later viewing or printing.

Requirements

COMPUTOR OS

Windows Operating System

NETWORK CONNECTIVITY

- Network connectivity for NeviFile computer
- Network connectivity for Nevisense [wireless is not supported]

The 10 Mbit Ethernet port on the Nevisense device shall only be used with Separation device in conformity with IEC/EN 60601-1:2005 §16.5.


Ordering information

SOFT WARE

142-0003 NeviFile kit including:

- NeviFile SW [provided on USB flash drive or as link]
- NeviFile Instruction for Use
- Network isolator device

Scibase AB
PO Box 3337, SE-103 67 Stockholm, Sweden
Phone: +46 8 410 620 00
E-mail: info@scibase.com
Order: order@scibase.com
www.scibase.com



Measure what can't be seen.